



Catlabs Pro 320

Test Date: 10/20/22
Curves: 2:45 min., 4 min., 5:45 min., 8 min., 11 min.
Film ISO: ISO 100
Development: D76 Kodak, 1+1, 68°F, 500
Agitation: Constant
Flare Density: none
Notes: Exposed on a DIY proofer/sensitometer




Catlabs Pro 320 / D76 1+1 / 20C




H&D curve for 5.75 min

Catlabs Pro 320

Test Date: 10/20/22
Curves: 2:45 min., 4 min., 5:45 min., 8 min., 11 min.
Film ISO: ISO 100
Development: D76 Kodak, 1+1, 68°F, 500
Agitation: Constant
Flare Density: 0.0200
Notes: Exposed on a DIY proofer/sensitometer



Catlabs Pro 320 / D76 1+1 / 20C



H&D curve for 5.75 min